

Disposable Circular Stapler

The Welfare Circular stapler is used in general surgery, thoracic surgery, surgical treatment of obesity and in colo-rectal surgery, to perform end-to-end, side-to-end, and side-to-side anastomoses.

Controlled tissue compression provides the flexibility to adapt to tissue depending on the clinical situation.

Adjustable Height staples accommodate compressed tissue thickness from 1.0mm to 2.5mm.



Disposable Circular Stapler

reference code	Number of Staples	Open Staple Height	Outer Diameter
WM-CS-23	20	4.8mm	23mm
WM-CS-26	24	5.3mm	26mm
WM-CS-29	26	5.3mm	29mm
WM-CS-32	30	5.3mm	32mm
WM-CS-34	32	5.3mm	34mm

